

Procedure for pancreatic lipase inhibitory activity assay (micelle method)

Written by Eisuke Kato on 2013/10/09

Modified on 2017/09/27

Modified on 2019/01/16

Reagents

1. Porcine pancreatic lipase (Sigma, L3126)
2. Glyceryl trioleate (Sigma, T7140-1G)
3. L- α -Phosphatidylcholine (Sigma, P4279, purity>99%)
4. Taurocholic Acid Sodium Salt (Nacalai tesque, 32729-74)
5. Cetilistat (Lipase inhibitor, TCI, C2745)
6. NEFA C-testwako (FUJIFILM Wako Pure Chemical Co., 279-75401)
7. Several solvents.

Preparation of reagent solution

A. Tris-HCl buffer

Prepare 13 mM Tris-HCl buffer, pH 8.0, with 150 mM NaCl and 3 mM CaCl₂)

B. Porcine pancreatic lipase solution

Dissolve 4.5 mg of Porcine pancreatic lipase in 30 mL of Tris-HCl buffer.

***Prepare fresh solution before the assay**

C. Glyceryl trioleate solution

Dissolve 1.0 g of Glyceryl trioleate in 6.25 mL of chloroform. Store under nitrogen atmosphere at -20 °C.

D. L- α -Phosphatidylcholine solution

Prepare 100 mg/mL ethanol solution. Store at -20 °C.

E. Stop reagent

Prepare 1 M aq. HCl.

F. Positive control (cetilistat)

Prepare 64 μ M solution in DMSO. Store at -20 °C. Dilute 2 times with water to prepare 32 μ M/50% aq. DMSO solution before the assay.

Procedure for pancreatic lipase inhibitory activity assay (micelle method)

Written by Eisuke Kato on 2013/10/09

Modified on 2017/09/27

Modified on 2019/01/16

Preparation of glyceryl trioleate emulsion

1. Add 200 μL of Glyceryl trioleate solution and 200 μL of L- α -Phosphatidylcholine solution in a glass tube and evaporate the solvent.
2. Add taurocholic Acid Sodium Salt (10 mg) and cold* Tris-HCl buffer (9 mL) in the tube.
*Ice cooled or stored in refrigerator.
3. Vortex the tube for 20 seconds and make a suspension
4. On ice, sonicate the suspension using ultrasonic homogenizer for 5 min and create a micelle solution (duty cycle: constant, output control: 3, time: 5 min, with Branson Sonifier 250)

*The prepared emulsion will be like the photo below.



Assay procedure

1. Add 200 μL micelle solution and 100 μL sample solution (in 50% DMSO aq.) to 1.5 mL micro tube, mix and pre-incubate for 5 min at 37 °C*
2. Add 100 μL of lipase solution, mix and start the reaction
3. After reacting for 30 min at 37 °C, add 40 μL of stop reagent. Mix well.
4. Add 600 μL of hexane
5. Mix rigorously for 60 sec. and leave settled until phase separation
6. Repeat step 5 for three times
7. Take 300 μL from hexane layer to the 0.6 mL micro tube
8. Evaporate the hexane and dissolve the residue in 100 μL of DMSO.
9. Measure the oleic acid in the DMSO solution by NEFA C-testwako kit.

Procedure for pancreatic lipase inhibitory activity assay (micelle method)

Written by Eisuke Kato on 2013/10/09

Modified on 2017/09/27

Modified on 2019/01/16

Required mixtures

	Micelle sol.	Sample sol.	Lipase	Buffer
Control	200 μ L	100 μ L	100 μ L	
Control Blank	200 μ L	100 μ L		100 μ L
Sample	200 μ L	100 μ L	100 μ L	
Sample Blank	200 μ L	100 μ L		100 μ L

Calculation

$$\text{Inhibitory activity (\%)} = \left(1 - \frac{A_{550\text{sample}} - A_{550\text{sample blank}}}{A_{550\text{control}} - A_{550\text{control blank}}} \right) \times 100$$